

## COLOUR RECIPES TO EMULATE GEL IN FRONT OF A TUNGSTEN SOURCE (3200 K)

Nando 502 Wash	Red	Green	Blue	Lime
	%	%	%	%
Open 3200 K Tungsten	85	21.4	11.8	100
G365 Warm Straw	100	7	5.6	85.6
L052 Light Lavender	61	0	19.4	81
L104 Deep Amber	100	0	0	86
L106 Primary Red	100	0	0	0
L115 Peacock Blue	0	100	16	14
L119 Dark Blue	0	68.4	100	0
L136 Pale Lavender	100	13.2	19.4	100
L161 Slate Blue	30	100	61.4	100
L174 Dark Steel Blue	100	100	42	84.2
L200 Double CTB	29.6	69	68.6	100
L201 Full CTB	100	85	34	100
L202 1/2 CTB	100	53.2	19.8	100
L203 1/4 CTB	97.6	38.6	15.2	100
L204 Full CTO	100	6.2	2.4	71
L205 1/2 CTO	90.4	5.2	4.8	100
L232 Correction White Flame Green	100	17.8	3.2	50
L236 HMI to tungsten	100	7.8	4	66.4
L728 Steel Green	100	87	18.8	100
L75 Evening Blue	16	86.4	80	37
L763 Wheat	70.8	5.6	7.4	100
R15 Deep straw	100	0	0	66.8
R21 Golden Amber	100	21.8	0	3.8
R3220 Double CTB	45.6	86	68.5	100
R371 Theatre booster 1	57.2	33.6	21.8	100
R385 Royal Blue	5.79	55.67	97.2	5.78
R76 Light Green Blue	0	100	21.6	0.4
R79 Bright Blue	2.8	69.6	41	0
R82 Surprise Blue	15.6	88.2	77.9	49.2
R85 Deep Blue	0	52.2	41	0